



# TRANSFORMATION USING CALCIUM CHLORIDE


---





# Objectives

- To understand the concept and procedure of calcium chloride transformation
- to successfully transform cells to constitutively express our gene of interest (GFP in a DH5alpha cell)





## BioBrick part(s)

- BBa1722008  
([http://partsregistry.org/wiki/index.php/Part:BBa\\_I13522](http://partsregistry.org/wiki/index.php/Part:BBa_I13522))
  - This is a constitutively expressed GFP which can be “turned off” with tetracycline
  - This part has been tested and proved to work

# Procedure


- Prepare competent cells
- Assess competency of cells
- Transform competent cells





## ***Materials***


- BioBrick Part: BBa\_I13522 pTet-GFP
- Single colony of *E. coli* cells: DH5alpha
- LB medium
- CaCl<sub>2</sub> solution, ice cold
- LB plates containing ampicillin
- Chilled 50-ml polypropylene tubes
- Beckman JS-5.2 rotor or equivalent
- 42°C water bath




# Prepare competent cells

- 1. **Inoculate** a single colony of *E. coli* cells into 50 ml LB medium. Grow overnight at 37°C with moderate shaking (250 rpm)
- 2. **Inoculate** 4 ml of the culture into 400 ml LB medium in a sterile 2-liter flask. Grow at 37°C, shaking at 250 rpm, to an OD 590 of 0.375.

- 3. **Aliquot** culture into eight 50-ml prechilled, sterile polypropylene tubes and leave the tubes on ice 5 to 10 min.
- 4. **Centrifuge** cells 7 min at  $1600 \times g$  (3000 rpm in JS-5.2),  $4^{\circ}\text{C}$ . Allow centrifuge to decelerate without brake.
- 5. Pour off supernatant and **resuspend** each pellet in 10 ml ice-cold  $\text{CaCl}_2$  solution.
- 6. **Centrifuge** cells 5 min at  $1100 \times g$  (2500 rpm),  $4^{\circ}\text{C}$ . Discard supernatant and **resuspend** each pellet in 10 ml ice-cold  $\text{CaCl}_2$  solution. Keep resuspended cells on ice for 30 min.
- 7. **Centrifuge** cells 5 min at  $1100 \times g$ ,  $4^{\circ}\text{C}$ . Discard supernatant and **resuspend** each pellet in 2 ml ice-cold  $\text{CaCl}_2$  solution.
- 8. **Dispense** cells into prechilled, sterile polypropylene tubes (250- $\mu\text{l}$  aliquots are convenient). **Freeze** immediately at  $-70^{\circ}\text{C}$ .






# Assess competency of cells


- 9. **Plate** appropriate aliquots (1, 10, and 25  $\mu\text{l}$ ) of the transformation culture on LB/ampicillin plates and incubate at  $37^{\circ}\text{C}$  overnight.
- 10. **Calculate** the number of transformant colonies per aliquot volume ( $\mu\text{l}$ )  $\times 10$



# Transform competent cells



- 11. 15 $\mu$ l of cell and 15 $\mu$ l of plasmaetes were taken in the tube and frozen in the ice bath for 30 minutes.
- 12. The LB broth solution was made by mixing 6.25g of LB broth powder in 250ml of DI water and autoclaving it for 15 minutes.
- 13. The tubes were placed in a water bath at 42°C for 90 seconds (heat shock step).
- 14. The tubes were cooled for 15 minutes in the ice bath.
- 15. 1ml of the prepared LB broth solution was added to the cells and shaken in the incubator for about an hour.
- Control was prepared by adding 1ml of the LB broth solution to the 15 $\mu$ l of the cells and shaking it in the incubator for about an hour.
- 15 $\mu$ l of the sample was transferred to the petridish and streaked with a sterilized glass rod.
- The petridishes were left in the incubator for about 12 hour for the growth.
- After 12 hours the petridishes were viewed under UV light to see the cells



# Outcome of experiment

- When the colonies were hit with UV rays from the UV lamp, in 2 out of the 4 petridishes, the cells exhibited a green fluorescent color
- When tetracycline was added to the medium, the cells no longer exhibited a green fluorescent color




BBa\_I13522 visualized  
under non-UV lightbox



BBa\_I13522 visualized  
under 254nm wavelength  
UV lightbox

# Skills learned in this lab exercise



- Transformation of Plasmid DNA into E.Coli cells using CaCl<sub>2</sub>, Heat Shock Transformation
  - Calcium ions help the uptake of the Plasmid DNA (making it competent)
  - Mixture of DNA and cells is heat shocked to allow entry of DNA into cells
  - Cells are grown to allow synthesis of proteins encoded in Plasmid DNA
- Lab technique