

β-Galactosidase Activity Assay (96-well plate format)

DAY 1:

1) Create solutions of LB + appropriate antibiotic (amp₅₀, amp₅₀+kan₃₅, tet₁₀)

e.g. LB+amp₅₀ solution

$$(4\text{columns} \times 4\text{rows}) + 1 = 17\text{mL LB}$$

$$17 \mu\text{L amp}_{50}$$

2) Put 1mL of appropriate solution in each well being used

Example: Shows 4 gates being tested by themselves and then with their 4 inputs, hence four columns of LB+amp₅₀ for the gates only, and four columns of LB+amp₅₀+kan₃₅ for gates+input and then in column 10 usually XL1-B is used for a blank so LB+tet₁₀.

LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+tet ₁₀		
LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+tet ₁₀		
LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+tet ₁₀		
LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+amp ₅₀ +kan ₃₅	LB+tet ₁₀		

3) Pick 4 colonies from plates containing appropriate plasmids and place 1 colony under each appropriate column in the four different rows.

4) Leave O/N to grow at 37°C and 320rpm.

DAY 2:

1) Create solutions of LB+appropriate antibiotics+IPTG+glucose/arabinose(2%)

e.g. LB+amp₅₀+IPTG+glucose

$$(4\text{columns} \times 4\text{rows}) + 1 = 17\text{mL LB}$$

$$17 \mu\text{L amp}_{50}$$

$$17 \mu\text{L IPTG}$$

$$1700 \mu\text{L glucose}(20\%)$$

2) Put 1mL of solutions into appropriate columns/rows of 96-well plate

LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+tet ₁₀ +IPTG+glucose		
LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+tet ₁₀ +IPTG+glucose		
LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+amp ₅₀ +kan ₃₅ +IPTG+glucose	LB+tet ₁₀ +IPTG+glucose		
LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+tet ₁₀ +IPTG+arabinose		
LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+tet ₁₀ +IPTG+arabinose		
LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+tet ₁₀ +IPTG+arabinose		
LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+amp ₅₀ +kan ₃₅ +IPTG+arabinose	LB+tet ₁₀ +IPTG+arabinose		

- 3) Inoculate 1mL with 100µL of each O/N (multi-channel pipettor)
 - e.g. row 1+5 same O/Ns
 - row 2+6 same O/Ns
 - row 3+7 same O/Ns
 - row 4+8 same O/Ns
- 4) Incubate at 37°C and 320rpm until A₆₀₀~0.3
- 5) Make ONPG 4mg/mL in 0.1M phosphate buffer (20µL/well)

Measuring Cell Concentration A₆₀₀:

- 1) Transfer 150µL of cells from each well into clear, flat bottom plate. (multi-channel pipettor)
- 2) Put LB in column 11 as blank
- 3) Read A₆₀₀

Lyse Cells:

- 1) Aliquot 1mL Z buffer + 20µL (0.1%) SDS + 40µL chloroform into a different, empty 96-well plate
 - *chloroform not in tray
- 2) Using multi-channel pipettor transfer 175µL of culture from other 96-well plate into this plate and pipette up and down 15 times to lyse cells
 - *Pipette from middle and dispense into bottom of the well where chloroform is
- 3) Leave for a while (5-10 minutes so chloroform settles to bottom)

LacZ Assay:

- 1) Transfer 200µL (middle to top) from lysed cells 96-well plate into clear 96-well plate.
- 2) As a blank, transfer 200µL of Z buffer into last column of clear plate
- 3) Add 20µL ONPG to all wells recording t₀
- 4) As colour develops add 50µL 1MNa₂CO₃, record t_f
- 5) Once all reactions are stopped use plate reader to record A₅₅₀(cell debris) and A₄₂₀(yellow product)

Calculations:

$$\text{Miller Units} = 1000 \times [(A_{420} - (1.75 \times A_{550}))] / (T \times V \times 4 \times A_{600})$$

A₄₂₀ and A₅₅₀ are read from the reaction mixture

A₆₀₀ reflects cell density in the washed cell suspension

T=time of the reaction in minutes

V=volume of culture used in the assay in mLs

Units give the change in A₄₂₀/min.mLL of cells/A₆₀₀