

Colony Screening PCR

Preparation – the day before

* White colonies are picked individually into 100ul LB/Carb in sterile 96 well plates (with a lid) and grown, standing, overnight at 37°C.

* The culture is used as template for PCR the next day (or may be held at 4°C prior to use).

* *Be sure to mix thoroughly prior to using the culture for csPCR.* If you don't, you may transfer too many cells and your PCR will not work.

Reaction Setup

Sign up for a PCR machine before you start-be considerate here if you are running lots of plates and take them out promptly.

Prepare a cocktail for all PCRs. Usually, n+1 reactions is sufficient, where n is your number of PCRs. However, where you may be setting up >10 reactions, you can use n+2 reactions for your cocktail, >50 reactions n+3, etc.

Reagent	Stock Concentration	Final Concentration	Vol]/20ul reaction
T7 primer	20uM	1uM	1.0 ul
SP6 primer	20uM	1uM	1.0 ul
dNTPs	2.5mM	200uM	1.5 ul
Buffer	10x	1x	2.0 ul
Taq	-	-	0.1 ul
H ₂ O	-	-	13.4 ul

Aliquot 19 ul per well in a 96-well plate. Add 1.0 ul of either overnight bacterial culture (feel free to use the 0.1-10 ul multichannel pipets for this) or positive control C2.

Reaction Conditions

Record the PCR machine you are using, and the positions in which your tubes are placed.

PCR Program: Method 14

Links programs 12, 13, 7, 8

94°C, 4 minutes.....Program 12

30 cycles:

94°C, 30 seconds

55°C, 30 seconds

72°C, 1 minute.....Program 13

72°C, 3 minutes.....Program 7

10°C, forever.....Program 8

*The appropriate extension time is 1 minute per 1kb of expected product size. Don't forget that you are adding on about 100bp to your PCR product due to flanking regions of the vector that are between the T7 and SP6 sites.