

Second device assembly		
	BBa_K119010:p.(TetR+luxR)+luxR/p.(luxR+ AHL)+cl*/p.lacZ+aiaA→Prk415 BBa_K11901:p.(TetR+luxR)+luxR/p.(luxR+ AHL)+cl*/p.lacZ*+aiaA→Prk415	
STEPS	DETAILS	PROBLEMS
Primer design	Our oligos were design specifically, so we could change GFP from Chiba biopart an put CI there Also, this oligos were bearing sites for restriction diqests. EcoRI-LuxR-Xbal -CI-HindIII -aiaa-PstI	
Extract DNA from the registry	Get alreay assembled biobricks we will use	But there were no DNA
Ask Chiba to send their biopart		We waited some time
Ask MIT for CI biopart		
PCR	Amplified bioparts and modified Chiba biopart Usinq specific oligonucleotides, getting 3 parts (LuxR and aiaa from Chiba) and CI from MIT	
Chek	Check our PCR product with electrophoresis on agarose gel	
Cloning into pJET	Blunting PCR product. pJET vector is already blunted Ligation:Usinq rapid ligation protocol, T4 DNA liqase and blunt-end ligation	
Transformation	Transformation: E. Coli - DH5a by heat-shock	
Check	Extraxt plasmid for selected colonies Restriction diqest of all colony extractions PCR with specific primers for our device Check the Restriction products and PCR on agarose gel	
Extraction	Plasmid extraction from confirmed colonies	
Purification	Doble restriction diqest of plasmids to extract our cloned device Gel band purification of restricted bioparts	
Assembly of the complete device	Ligation of the three individual devices Sticky end liqations	We just be able to liqate two of three parts One oliqo got bound somewhere in aiaa biopart.
Cloning into pBBR1MCS-5	Restriction Digest of pRK415 vector Ligation: pRK415 + first bioparts Ligation: prk415 + liqation of two bioparts	We had a lot of troubles trying to clone our bioparts in PRK
Transformation	Transformation: E. Coli - DH5a by heat-shock	
Extraction	Extract plasmid for selected colonies (XqaI control)	Very few white colonies.
Check	Restriction diqest of all colony extractions PCR with specific primers for our device Check the Restriction products and PCR on aqarose qel	Many false positives with XqaI control
Alternatives (Not realized)	Ligate all the bioparts in one, using a rapid ligation protocol, keeping the sites for pJet oligos Hybrid plasmid: a plasmid made by the joint of two different plasmids, pJet and PRK with DNA ligation using shrimp alkaline phosphatase	Not done
Transformation of the final receptor cell	Transformation of E. Coli Yoh- with the assembled device	Not done
Check	Extraxt plasmid for selected colonies Restriction diqest of all colony extractions PCR with specific primers for our device Check the Restriction products and PCR on aqarose qel Add suffix and prefix Clonate the device in a standard plasmid Get sequence of our final construction	Not done