

Wrapup 2008-06-24 1 of 2

Priority

General

Plasmids

pLV-TRE-Ngn1-EYFP-2A-mKate-Ubc-Bla

[1] O.D. 192 ng/ul for Ngn1-EYFP-2A-mKate; ~300ng/ul for Ngn1-EYFP

[2] Insert digest with SfiI

[2] PCR purify insert digest.

[2] Measure O.D.

[1] Digest pLV-TRE-Sox17-Ubc-Bla with SfiI [double-check that]

[1] Run vector digest on gel.

[1] Redo digest – 10ug (or 20ug if we have it). [David will re-digest tonight].

[1] Run on 1% agarose gel at a higher voltage (110V) for 2 hours on a longer gel (no middle combs). [try to get rid of u-shaped bands – ask VectorNTI]

[1] Gel extract. (if not good separation, take the bottom part of the band.)

[1] Measure O.D.

[2] CIP vector digest

[2] PCR purify vector digest CIP

[2] Measure O.D.

[2] Ligate

[3] Transform

pLV-TRE-Sox17-Ubc-Bla

[3] Grow two colonies overnight [Evan and Navin]

[3] Midiprep. – O.D. 250 and 270ng/ul

[3] Restriction map with BamHI and MluI; expect bands at 7498 and 3461.

pLV-Ubc-rtTA-2A-Bla

[3] Sunday Evan and Lena will pick 40 minipreps (10 from each plate) for growth.

[3] Extract DNA

[3] O.D. 30-40ng/ul

[2] Design restriction map.

[2 and 3] Restriction map – cut with MluI. Do as many as you can as soon as you can.

[3] If you find a good one – transform today to maxiprep tomorrow

[3] On the three decent-looking samples [7-2, 7-8, 7-6], transform, midiprep, sequence, restriction map again.

p148

[3] O.D. [~500 ng/uL]

p149

[3] Maxiprep [1330-1360 ng/ul] (Two O.D. tries)

pLV-TRE-PKMz-Ubc-Hyg

[1] Design plasmid

Wrapup 2008-06-24 2 of 2

[1] Buy PKMz

pLV-pPKMz-EGFP-Ubc-Bla

[1] Design Plasmid

[1] Buy the promoter (PKMz)

pLV-TRE-Mash1-EYFP-2A-mKate-Ubc-Bla **(in same plasmid as Ngn1?)**

[1] Design plasmid [do this one first]

pLV-TRE-Mash1-Cerulean

pLV-TRE-Ngn1-EYFP

pLV-NeuronalPromoter-Neuronal specific CFR

??

Lentivirus

Lenti: pLV-TRE-Ngn1-EYFP-2A-mKate-Ubc-Bla

Experiments

Experiment: Ngn1-EYFP-2A-mKate

Infect cells

Add Dox

Observe differentiation

Make Yellow/Red artificial brains

Plan ahead: Camera sensitivity, Dyes, Cil found a dye which could help us see action potentials.

Buy some neurons, Dyes, potassium and valinomycin, glutamate, GABA, dopamine, ACh, special media

Start planning microfluidics and also the optical twizzlers

Bio nanoforce