

# Wrapup 2008-07-09 1 of 3

## Possible Plasmids:

pLV-TRE-Ngn1-EYFP-2A-mKate-Ubc-Bla:

A6 Midiprepped, PCR purified, sequence tonight; 68 ng/ul

B4 Midiprepped, PCR purified, sequence tonight; 45 ng/ul

A22 Grow for midiprep: **Ngn1 sequencing successful! ;**

A23 Midiprepped, PCR purify tomorrow; 62 ng/ul

A25 (sequenced) Excess DNA????

B14

B18

## Priority

## General

## Debugging

## Plasmids

pLV-TRE-Ngn1-EYFP-2A-mKate-Ubc-Bla

[3] Minipreps – 25 from each plate (50 total)

[2] O.D. – 4 – 108ng/ul; mostly 20-60ng/ul

[2] Finish Restriction digests with AscI and BsmBI

[2] Run restriction digest on gel for last 18. Restricted with Bsu36I; expect bands at ~3000 and ~9000bp.

[1] Sequence new Ngn1 – Designed and ordered

[1] Sequence insert – A22, A25, [David is taking it down tonight]

[3] Midiprep (with smaller elution volume) A6 (54ng/ul), B4 (17ng/ul)

[2] PCR purify maximum volume of A6 (68 ng/ul), B4 (45 ng/ul) that fits into a single tube [100 or 120ul?], elute with 40ul

[2] PCR purify the DNA with higher OD (A23)

Molly and Eric will find a good concentrator to use and will be very kind to the head of said lab so that we can use it the rest of the summer.

[1] Sequence A6, B4, rtTA [Andrew]

[1] A22 sequence of Ngn1 successful by BLAST; ALIGNX verification in Vector NTI

[3] Maxiprep A22

[1] Sequence A22 Maxiprep properly

[1] Design/order correct sequencing primer (instead of UBC forward, order a reverse primer)

[3] Grow and midiprep A23 (Transform, Grow, Midiprep)

[3] Midiprep B14, B18 - transform

[2] Restriction digest A6, B4 with Bsu36I (the ones that looked promising but hadn't been digested with Bsu36I yet)

pLV-TRE-Sox17-Ubc-Bla

[1] Sequencing (use existing primers)

pLV-Ubc-rtTA-2A-Bla

[1] Sequence (wait for primers)

[3] Restriction map with MluI – expect bands at 6386 and 3933 – gel is running.

[3] Restriction digest and map again with remaining minipreps

# Wrapup 2008-07-09 2 of 3

[1] Find a single-cut enzyme - PacI

[1] Digest

[1] restriction map

Bad sequencing results.

[2] PCR purify 20ug or 120ul, whichever is larger

[1] Design new sequencing primers

[1] Sequence the PCR purified DNA

pFUGW

[3] Make more pFUGW – transformed

[3] Maxiprep – O.D. = 250ng/ul

p148

[3] O.D. [~500 ng/uL]

p149

[3] Maxiprep [1330-1360 ng/ul] (Two O.D. tries)

pLV-TRE-PKMz-Ubc-Hyg

[Cil] Request gene from lab X

[1] Design plasmid

[1] Buy PKMz

pLV-pPKMz-EGFP-Ubc-Bla

[Cil and Team #1] PCR PKMz gene and promoter out of chromosome

[1] Design Plasmid

[1] Buy the promoter (PKMz)

pLV-Hef1a-LacO-Mash1-Cerulean-Ubc-Hyg

[1] Show design to Cil

[1] Order primers

[1] PCR SOEing to make the plasmid – add Ubc-Hyg to PCR SOEing

Design SOEing primers for Mash1-Cerulean-Ubc-Hyg

[1] The backbone's going to be Hef1a/LacO

pLV-TRE-Ngn1-EYFP-Ubc-Bla

[1] Gradient PCR EYFP by itself from the original plasmid [Eric and Caroline Jr.]

[1] Run gel

[3] Transform and midiprep pEYFP

[1] Design digest (50-250 bp on either side of EYFP) and digest pEYFP with (5') HindIII [37C] or BsaXI and (3') SpeI or BsmBI [55C]

Verify result of enzyme check.

[1] PCR EYFP

[1] PCR Ngn1

[1] SOE Ngn1-EYFP.

[1] Run gel

[1] Extract

pLV-NeuronalPromoter-Neuronal specific CFR

# Wrapup 2008-07-09 3 of 3

??

## Vector NTI Stuff

What we have in the database / load database

Look over primer design for Ngn1-EYFP

Put primer design manual in public folder

Design PCR SOEing primers for pLV-Hef1a-LacO-Mash1-Cerulean-Ubc-Hyg

## Lentivirus

Lenti: pLV-TRE-Ngn1-EYFP-2A-mKate-Ubc-Bla

## Experiments

Experiment: Ngn1-EYFP-2A-mKate

Infect cells

Add Dox

Observe differentiation

Make Yellow/Red artificial brains

Plan ahead: Camera sensitivity, Dyes, Cil found a dye which could help us see action potentials.

Buy some neurons, Dyes, potassium and valinomycin, glutamate, GABA, dopamine, ACh, special media

Start planning microfluidics and also the optical twizzlers

Bio nanoforce